



## Single-Cell Perturb and TAP Sequencing

### Product overview

**Perturb-seq** is a high-throughput pooled CRISPR screen that combines CRISPR technology for precise gene manipulations with single-cell RNA sequencing (10x Genomics) allowing for whole transcriptome profiling of perturbation effects.

**TAP-seq** (TARgeted Perturb-seq) has a similar approach but focuses on a subset of genes, selected by the user, using gene-targeting primer pools. TAP-seq is cost effective and is more sensitive for lowly expressed genes.

Sequencing will be performed on the NovaSeq X Plus platform (Illumina) with tailored sequencing settings based on the project design and the user's needs.

### What is included

- GEM generation
- Single-cell library preparation
- Quality control of cDNA and library
- Library pooling and sequencing on Illumina NovaSeq X Plus
- FASTQ files
- Standard data processing, QC, and basic bioinformatics analysis (optional).

### Input and quality requirements

#### Isolation of cells or nuclei

Preparation and thawing of perturbed cells as well as cell counting is the responsibility of the user. Contact [SCOP](#) to consult.

#### Quality

The quality of the cells will reflect the downstream process – poor quality leads to poor data.

- > 90% cell viability (a cell viability between 70-80% can be accepted but data quality may be compromised)
- Intact cell membranes
- No to minimal debris and cells aggregates.

SCOP strongly recommends performing a pilot experiment to assess both sample and data quality.

#### Sample drop-off

Perturbed cells should be delivered at a specific concentration which depends on the assay type and the desired number of cells. Contact [SCOP](#) for more information.

### Data deliverables

SCOP's preprocessing workflow includes:

- FASTQ files.
- Preprocessed data stored in Seurat objects.
- Output from Sceptre.