



## Spatial Transcriptomics

### Product Overview

The **10x Genomics Xenium Analyzer** is a high throughput and flexible In Situ Platform which enables high-plex in situ analysis at subcellular resolution on fresh frozen (FF) and formalin fixed paraffin embedded (FFPE) tissue. Xenium allows the following workflows:

#### Xenium In Situ Gene Expression (v1):

Targets up to 480 gene targets per slide and has a variety of probe panel options which makes the workflow very flexible; these include: pre-designed panels, custom add-ons, standalone custom or advanced custom panels.

#### Xenium In Situ Gene and Protein Expression:

An extension to the v1 workflow targeting both RNA and protein simultaneously. Protein is detected through one of six available antibody subpanels and is restricted to FFPE tissue sections.

#### Xenium Prime 5k Gene Expression:

Targets 5000 genes per slide with the option of adding a custom add-on probe panel (up to 100 custom gene targets). The Prime workflow is limited to two species: human and mouse.

The Xenium Analyzer utilizes DAPI-staining to perform cell segmentation by nuclear expansion. To optimize cell segmentation, users can either include the Xenium Cell Segmentation Add-on Kit (during the Xenium run), and/or perform manual staining (e.g. IF-staining) after the Xenium run. Please find a list of tissues validated using the Xenium Cell Segmentation kit [here](#). If performing manual segmentation or optimization of segmentation is required, Cellpose, a free Deep learning-based image segmentation tool, can be used.

### What is Included

- Tissue detachment test (optional)

- Xenium workflows (V1, V1&protein, or Prime 5K)
- Xenium run
- Data delivery
- Post-stain image alignment upon request.

### Input Requirements

Xenium slides with high-quality tissue sections:

- Fresh frozen (FF)
- Formalin fixed paraffin embedded (FFPE)
- Tumor microarrays (TMAs).

### Gene panels

Choose a panel that is suited for your research project: 10x Genomics pre-designed panels, custom add-ons, standalone custom, advanced custom panels (e.g. other species).

If add-on or custom probe panels are required, design and order probes through 10x Genomics (due to IP rights). Remember to send the probe *.son* file to SCOP.

### Data Deliverables

- Xenium Ranger output files
- Re-segmented with a custom Cellpose model upon request.
- Segmentation refinement using Baysor
- Combined clustering and cell annotation across slides
- Generation of an annotated data object.

### Disclaimer

SCOP does not offer the following services: preservation/embedment of tissue, tissue sectioning and/or mounting (contact RMPP), post-stain (e.g. IF-staining, H&E-staining), building cell segmentation models, custom data analysis.